	control pool CV intrapla de (median)	control pool CV interpla de	median LLOD (1:10)	ULOD (1:10)	samples in LOD (n)	samples in LOD %	LLOQ (1:5)	ULOQ (1:5)	samples in LOQ (n)	samples in LOQ (%)	out of range (n)	out of range (%)
Osteopontin	26,4	27,1	422	500000	9196	99,9	8000	528979	9134	99,2	8	0,1

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
units for all

ULOD: 0 std. above S001

LLOQ: lowest standard which has CV <25% and %recovery 75-125%

ULOQ: S001, if it has CV <25% and %recovery 75-125%

MSD method description **2-plex Follistatin and MPO**

The concentrations of Follistatin and MPO were measured on a MESO Quickplex SQ 120 instrument (Mesocore diagnostics (MSD)), using R-plex assay kits (F21E5-3/-8, F214E-3/-8) according to the manufacturer's instructions. The assay is based on analyte specific antibody pairs, with a primary antibody that captures the analyte to one of up to 10 different spots in a 96well and a secondary antibody that generates a light signal proportional to the analyte amount upon electrical stimulation. Light signals are converted to concentrations using a standard curve, which is generated for each individual plate.

9204 samples were measured on 118 individual plates, each plate containing 78 samples in single wells, as well as 8 standards and a control pool in duplicate. The control pool was used to determine the average assay variation, both intra and inter-plate.

Samples underwent one freeze thaw cycle to be aliquoted onto 96 well plates before they were thawed in the 96 well plates for the MSD analysis. Thawed samples were spun down and diluted 1:1 with assay diluent for the analysis. All pipetting steps were performed with the BRAVO pipetting platform.

The average assay variation within and between plates was 2,7% and 7,9% for Follistatin and 3,9% and 15,8% for MPO. 99,95% and 99,98% of samples were measured in detection range of the assay for Follistatin and MPO, respectively.

Detection range was calculated for every plate based on the highest and lowest standard.

Units is given as pg/ml.

MSD method description **Adiponectin**

The concentration of Adiponectin was measured on a MESO Quickplex SQ 120 instrument (Mesocale diagnostics (MSD)), using a R-plex assay kit (# K151YTR, MSD) according to the manufacturer's instructions. The assay is based on analyte specific antibody pairs, with a primary antibody that captures the analyte to one of up to 10 different spots in a 96well and a secondary antibody that generates a light signal proportional to the analyte amount upon electrical stimulation. Light signals are converted to concentrations using a standard curve, which is generated for each individual plate.

9204 samples were measured on 118 individual plates, each plate containing 78 samples in single wells, as well as 8 standards and a control pool in duplicate. The control pool was used to determine the average assay variation, both intra and inter-plate.

Samples underwent one freeze thaw cycle to be aliquoted onto 96 well plates before they were thawed in the 96 well plates for the MSD analysis. Thawed samples were spun down and diluted 1:10.000 for the analysis. All pipetting steps were performed with the BRAVO pipetting platform, sample plates were incubated on an orbital shaker between the dilution steps.

The average assay variation within and between plates was 5,1% and 10,1%, respectively. 99,9 % of samples were measured in quantification range of the assay.

Quantification range was defined by the concentration of standards that had a CV <25% and a recovery of 75-125%.

Units is given as pg/ml.

MSD method **R-plex Osteopontin**

The concentration of Osteopontin was measured on a MESO Quickplex SQ 120 instrument (Mesocale discovery (MSD)), using a R-plex antibody set (K151YMR-2) and 96-well Smallspot Streptavidin plates (L45SA-1) according to the manufacturer's instructions. The assay is based on analyte specific antibody pairs, with a capture antibody that binds to a spot in a 96-well and a detection antibody that enables light generation, proportional to the amount of bound analyte. Light signals are converted to concentrations using a standard curve, which is generated for each individual analyte on each individual plate.

9204 samples were measured on 118 individual plates and each plate contained 78 samples in single wells, as well as 8 standards and a control pool in duplicate. The control pool was used to determine assay variation, both intra and inter-plate.

Samples underwent one freeze thaw cycle to be aliquoted onto 96 well plates before they were thawed in the 96 well plates for the MSD analysis. Thawed samples were spun down and diluted 1:10. All pipetting steps were conducted with the BRAVO pipetting platform.

Samples from x plates were fitted with the standard curve from another plate.

The average assay variation within and between plates was: x % of samples were measured in detection range of the assay.

Detection range was calculated for every plate based on the highest and lowest standard and quantification range was given by the supplier (V-plex) / was defined as the standards that had a CV <25 and a recovery of 75-125%.

	control pool CV intraplade (median (avera control pool CV interplade	median LLOD	ULOD	samples in LOD (n)	samples in LOD %	LLOQ	ULOQ	samples in LOQ (n)	samples in LOQ (%)	out of range (n)	out of range (%)	
Osteocalcin-1 (alle)	4,6 (8,9)	11,3	74,9	200000,0	9203,0	99,99	220,0	190370,0	9203,0	99,99	1,0	0,0
Osteocalcin-1 (coat 1h)	4,2 (4,3)	10,0	69,9	200000,0	9203,0	99,99	220,0	190370,0	6084,0	100,00	0,0	0,0
Osteocalcin-1 (coat 24h eller batch)	7,8 (17,9)	13,6	90,4	200000,0	9203,0	99,99	726,0	13478,0	529,0	16,95	0,0	0,0

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



LLOD: 2.5 std. above S008 (technical blank)

ULOQ: 0 std. above S001

LLOQ: lowest standard which has CV <25% and %recovery 75-125%

ULOQ: S001, if it has CV <25% and %recovery 75-125%

General information on assay/data quality:


Definitions:

LLOD: 2.5 std. above S008 (technical blank)

ULOD: 0 std. above S001

LLOQ: lowest standard which has CV <25% and %recovery 75-125%

ULOQ: S001, if it has CV <25% and %recovery 75-125%

Result Excel file contains column "RANGE" indicating whether a sample lies outside detection range "0", within detection range "1" ($LLOD < \text{sample} < LLOQ$) or within quantification range "2" ($LLOQ < \text{sample} < ULOQ$).

Overview Vplex panel IL6 and TNF α in DD2 cohort


TNF α and IL6 were detected with a Vplex assay, which stands for validated assay. This type of assay undergoes extensive testing and quality control by the manufacturer, which results in higher reproducibility, broader dynamic range, higher sensitivity, precision, recovery and linearity, compared to the other available assay types.

The TNF α assay was very stable and almost all samples were within the detection range, about 90% in the quantification range.


IL6 assay showed more variability, standards in the lower end of the curve had slightly elevated recovery values, which means that the measured signal/concentration was above the calculated standard curve. Due to this higher variation in the lower end (and of course the biologically given distribution of IL6 concentrations), fewer samples are within detection range and even fewer in quantification range. With very stringent criteria the lower limit of the quantification range should be between standard 5 and 6, since the criterion for lower quantification limit is only met in 67% of the measured 118 plates.

The reason for this higher standard recovery% seems to arise from slightly different standard preparation/dilution, apparent as overall lower standard signal in the plates where recovery is elevated. This is supported by TNF α standard signal values, but the assay seems to be more robust since recovery % is not affected here.

Standard recovery

CV on standard measurements

IL6

Detection limits (pg/ml)

- Median calculated LLOD 0,151 (min 0,10 max 0,35)
- ULOD 706

Quantification limits

Kit

- LLOQ_{KIT} **0,633**
- ULOQ_{KIT} **488**

Variation on internal control measurements

- Interplate CV on control pool concentration: 12,34
- Median of Intraplate CV on control pool concentration: 6,08
- Control recovery on the individual plates, relative to median control pool concentration (80,4-148,4%)

Number of samples within range

within detection range (according to limits on the respective plate)
9060

within quantification limits (*Kit*) (0,633< samples<488)
8015

% detected (I intervallet LLOD-ULOD) 98,5
% quantifiable (I intervallet LLOQ-ULOQ KIT) 87,1
% not detected (under LLOD) 1,5

Distribution of calculated concentration values (pg/ml):

MAX	191,9
MIN	0,04
AVERAGE	1,8
MEDIAN	1,2

TNFalpha

Detection limits (pg/ml)

Median calculated LLOD **0,032** (min 0,02 max 0,10)

ULOD **358**

Quantification limits (pg/ml)

KIT

Kit LLOQ 0,690

Kit ULOQ 248

Variation on internal control measurements

Interplate CV on control pool concentration: 14,4%

Median of Intraplate CV on control pool concentration: 4,98%

Control recovery on the individual plates, relative to median control pool concentration (72,8-135,1%)

Number of samples within range

within detection range according to limits on the respective plate

9201 (2 below curve fit, including the empty sample well, 1 below detection range)

within quantification limits (KIT) (0,69< samples<248)

8229

- % detected (I intervallet LLOD-ULOD) 99,99%
- % quantifiable (I intervallet LLOQ-ULOQ_{KIT}) 89,4%
- % not detected (under LLOD) 0,01%

Distribution of calculated concentration values:

MAX	86,4
MIN	0,04
AVERAGE	1,1
MEDIAN	0,98

MSD method **R-plex Osteopontin**

The concentration of Osteopontin was measured on a MESO Quickplex SQ 120 instrument (Mesocale discovery (MSD)), using a R-plex antibody set (K151YMR-2) and 96-well Smallspot Streptavidin plates (L45SA-1) according to the manufacturer's instructions. The assay is based on analyte specific antibody pairs, with a capture antibody that binds to a spot in a 96-well and a detection antibody that enables light generation, proportional to the amount of bound analyte. Light signals are converted to concentrations using a standard curve, which is generated for each individual analyte on each individual plate.

9204 samples were measured on 118 individual plates and each plate contained 78 samples in single wells, as well as 8 standards and a control pool in duplicate. The control pool was used to determine assay variation, both intra and inter-plate.

Samples underwent one freeze thaw cycle to be aliquoted onto 96 well plates before they were thawed in the 96 well plates for the MSD analysis. Thawed samples were spun down and diluted 1:10. All pipetting steps were conducted with the BRAVO pipetting platform.

Samples from x plates were fitted with the standard curve from another plate.

The average assay variation within and between plates was: x % of samples were measured in detection range of the assay.

Detection range was calculated for every plate based on the highest and lowest standard and quantification range was given by the supplier (V-plex) / was defined as the standards that had a CV <25 and a recovery of 75-125%.

(this was page 2 in the document from adiponectin)